



GUARANTEED
3
YEARS

AEROSOL CON COMPRESSORE A PISTONE PER USO PROFESSIONALE PISTON COMPRESSOR NEBULIZER FOR PROFESSIONAL USE

DATI TECNICI COMPRESSORE

COMPRESSORE	A PISTONE
PRESSIONE MASSIMA	280 kPa
PRESSIONE D'ESERCIZIO	0 ÷ 110 kPa
PORTATA D'ARIA	15 lt/min
CICLO DI UTILIZZO	CONTINUO
RUMOROSITÀ	55 dBA
DIMENSIONI	19,6 x 35,7 x 18,5H cm
PESO	2,5 kg
CLASSE DI RISCHIO	Ila

COMPRESSOR TECHNICAL DATA

COMPRESSOR	PISTON COMPRESSOR
MAX PRESSURE	280 kPa
OPERATING PRESSURE	0 ÷ 110 kPa
AIR FLOW	15 lt/min
OPERATING CYCLE	CONTINUOUS
NOISE LEVEL	55 dBA
DIMENSION	19,6 x 35,7 x 18,5H cm
WEIGHT	2,5 kg
MDD CLASS	Ila



MASCHERA ADULTO - ADULT MASK



TUBO ARIA - AIR TUBE



AMPOLLA FASTERJET - FASTERJET NEBULIZER



FILTRO POREX - FILTER POREX



MASCHERA PEDIATRICA - CHILD MASK



NASALE - NOSEPIECE



BOCCAGLIO CON VALVOLA
MOUTHPIECE WITH VALVE



DOCCIA NASALE - NASAL IRRIGATION

CE0434

(dir. 93/42/EEC)

PCS. **2**

MASTER CARTON

ACCESSORI

- AMPOLLA FASTERJET CON SISTEMA VALVOLARE (POLIPROPILENE)
- DOCCIA NASALE (POLIPROPILENE)
- NASALE (POLIPROPILENE)
- BOCCAGLIO (POLIPROPILENE) CON VALVOLA (SILICONE)
- MASCHERA ADULTO (PVC DOP FREE)
- MASCHERA PEDIATRICA (PVC DOP FREE)
- TUBO ARIA (PVC DOP FREE)
- FILTRO DI RICAMBIO
- ELASTICI

ACCESSORIES

- FASTERJET NEBULIZER WITH VALVE SYSTEM (POLYPROPYLENE)
- NASAL IRRIGATION (POLYPROPYLENE)
- NOSEPIECE (POLYPROPYLENE)
- MOUTHPIECE (POLYPROPYLENE) WITH VALVE (SILICONE)
- ADULT MASK (PVC DEPH FREE)
- CHILD MASK (PVC DEPH FREE)
- AIR TUBE (PVC DEPH FREE)
- SPARE FILTER
- ELASTIC

DATI TECNICI AMPOLLA

CAPACITÀ	16 ml
FRAZIONE RESPIRABILE	< 5 µm 80%
MMAD (Mass Median Aerodynamic Diameter)	2,7 µm
NEB RATE	0,60 ml/min
VOLUME RESIDUO	0,7 ml

NEBULIZER TECHNICAL DATA

CAPACITY	16 ml
RESPIRABLE VOLUME	< 5 µm 80%
MMAD (Mass Median Aerodynamic Diameter)	2,7 µm
NEB RATE	0,60 ml/min
RESIDUAL VOLUME	0,7 ml